PharmaJet®

Stratis®

EN-Needle-free Injection System

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CE
2797

Symbols Glossary

(Note: All symbols are derived from **ISO 15223-1, Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied.**)

Symbol	Title of Symbol	Description of Symbol	Ref. No.
	Manufacturer	Indicates the medical device manufacturer.	5.1.1
	Authorized representative in the European Community/ European Union	Indicates the authorized representative in the European Community/ European Union.	5.1.2
	Use-by date	Indicates the date after which the medical device is not to be used.	5.1.4
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	5.1.5
	Catalogue number	Indicates the manufacturer's catalogue number so the batch or lot can be identified.	5.1.6
	Date of manufacture / Country of manufacture	Indicates the date when the medical device was manufactured / to identify the country of manufacture of products.	5.1.3/ 5.1.11
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.	5.2.3
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.	5.2.8
	Single sterile barrier system	Indicates a single sterile barrier system.	5.2.11
	Single sterile barrier system with protective packaging outside	Indicates a single sterile barrier system with protective packaging outside	5.2.14
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	5.3.7
	Do not re-use	Indicates a medical device that is intended for one single use only.	5.4.2
	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use.	5.4.3
	Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.	5.4.4
	Non-pyrogenic	Indicates a medical device that is non-pyrogenic.	5.6.3
	Medical Device	Indicates the item is a medical device.	5.7.7
	Unique device identifier	Indicates a carrier that contains unique device identifier information.	5.7.10
	Prescription Only	Caution: US Federal law restricts this device to sale by or on the order of a physician.	N/A

PharmaJet®

Stratis®

English

System Components

Injector

Claimed Lifetime: 20,000 cycles

Reset Station

Syringe

STERILE EO

Filling Adapter

STERILE EO

CAUTION: Federal law restricts this device to sale by or on the order of a physician.

Indications for Use

The PharmaJet **Stratis**® Needle-free Injection System is intended to deliver various medications and vaccines either intramuscularly or subcutaneously by means of a narrow, high velocity fluid jet, which penetrates the skin and delivers the medicine or vaccine to the body. Healthcare providers who routinely administer injections may use the PharmaJet Needle-free Injection System. It may be used for adults and children. It can also be used by patients authorized by their physicians to self-inject, or have other individuals administer injections of prescribed medication.

Clinical Notes

The PharmaJet Stratis NFIS can be used with medications intended for injection into the intramuscular or subcutaneous tissue. Types of medications that may be used with the Stratis NFIS include anti-infectives (e.g., vaccines), antineoplastic and immunomodulating agents, hormonal preparations, and other therapeutic products intended for 0.5mL intramuscular or subcutaneous injection. It is the responsibility of the prescribing physician to ensure medicines are used within their licensed indications and methods of administration.

The following local responses may occur with Needle or Jet Injection Technology.

- a) Some medications may produce a burning, aching, or stinging sensation. Check the medication information for a description of local reactions.
- b) Following the injection, a drop of blood or fluid may appear. Blood and/or fluid loss can be minimized by immediately applying pressure to the injection site for approximately 60 seconds with a cotton ball or gauze. Apply adhesive bandage if needed or desired.
- c) To ensure delivered dose to target tissue, refer to **Section 4) Give Injection** for proper injection technique. For a partial or incomplete injection, refer to the medication or vaccine product labeling.
- d) In the event of changes in performance of the device, DO NOT USE, contact PharmaJet.

Cleaning, Disinfecting, and Storage

CAUTION: DO NOT autoclave.


CLEANING & DISINFECTING: Injector and Reset Station may be cleaned by wiping excess debris from the exposed surfaces using 70% Isopropyl Alcohol.

STORAGE: To store Injector in Reset Station, insert Injector and lift upward on the two sliding arms while closing the lid. This action allows Injector to be stored in the released state.

NOTE: The container and Injector are resistant to short term duration exposure to water.

DISPOSAL: Dispose of used disposables in accordance with local requirements.

Storage and Transport Temperature:


Reorder Information

<u>Item</u>	<u>Catalog #</u>
Injector	C501-10
Reset Station	C501-20
Syringes – 100 Bulk units	C501-34
Universal Filling/Vial Adapters – 100 Single units	C501-43
13mm Filling/Vial Adapters – 100 Single units	C501-46

CAUTION: Use only PharmaJet components with the Stratis® Needle-free Injection System.

Warranty Information

The Injector and Reset Station carry a replacement warranty against defects in materials and workmanship for two years from the date of purchase.

Contact Information

PharmaJet®

US Office:

PharmaJet, Inc.
400 Corporate Circle, Suite N
Golden, Colorado 80401
USA

Tel: +1 888 900 4321
Tel Outside the U.S.: +1 303 526 4278

EC	REP
----	-----

EU Authorized Representative:

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Suitable Injection Sites

WARNING: DO NOT position the Syringe directly over a blood vessel. Although injection into the vessel is highly unlikely, injection over a vessel may produce an adverse reaction including, but not limited to, bleeding, bruising, swelling, or tenderness.

Select a site that is free from irritation, with a low likelihood of infection, and where the skin is intact and free from abrasions, bruises, etc.

Needle-free injection sites are identical to those used for conventional needle and syringe; the deltoid, triceps and thigh areas are most often used as the delivery sites for vaccine. Other types of injectant may have specific sites for injection. Follow package directions from the manufacturer of injectant.

Warnings and Cautions

1) Prepare Injector

 **CAUTION:** Repeated dry firing will damage Injector. Dry firing is activating the Injector with either an empty Syringe or without a Syringe.

 **CAUTION:** The Injector may be damaged if not fully nested in the Reset Station when closing the lid.

2) Fill Syringe

 **WARNING:** The tip of the Filling/Vial Adapter is sharp. Use caution when handling.

 **WARNING:** Syringes auto-disable during use. Attempting to use a Syringe more than one time may result in contamination or device malfunction.

 **WARNING:** Vial integrity or sterility is not ensured when used with the Filling/Vial Adapter. **DO NOT** reuse the Filling/Vial Adapter in another vial. Doing so may contaminate contents of new vial.

 **CAUTION:** The protective Syringe Cap is designed to protect against touch contamination only, and does not provide a barrier against other forms of contamination.

 **CAUTION:** Visually inspect sterile package for damage. **DO NOT** use if torn or damaged. **DO NOT** re-sterilize.

 **CAUTION:** **DO NOT** attempt to use a partially filled Syringe; damage may occur.

3) Load Injector

 **WARNING:** **DO NOT** direct the Syringe toward the face or nearby skin when loading the Syringe into an Injector, preparing to give an injection, or ejecting a Syringe after use. Injury may occur.

 **CAUTION:** Partial or incorrect insertion of a Syringe into an Injector may damage the Syringe and/or Injector.

 **CAUTION:** **DO NOT** depress the Multi-function Button if the Polished Core Ring is not visible as shown in figure 10. This may cause inadvertent activation.

4) Give Injection

WARNING: Improper injection site location or technique may result in injury.

CAUTION: Patients receiving anticoagulants or persons with hemophilia or other coagulation disorders may have more post-injection bleeding and/or bruising.

CAUTION: If the Injector fails to reset or if the retractable core does not move freely, **DO NOT USE**. Continued use may result in an incomplete injection or failure to inject.


CAUTION: **DO NOT** deliver injection with the Injector at an angle other than 90 degrees to the patient.

Component Description


Injector

1. Multi-function Button
 - Safety
 - Activation
2. Polished Core Ring
3. Injector Hammer
4. Syringe Ejection Button


Reset Station

1. Reset Station Lid
2. Reset Station Base

Syringe

1. Syringe Cap
2. Syringe Barrel
3. Syringe Plunger
4. Fill Line

Filling/Vial Adapter

1. Syringe Interface
2. Filling/Vial Adapter Spike

Workflow

1) Prepare Injector

NOTE: Do not dry fire the Injector, damage will occur and the device may become unusable.

- a) Ensure the Reset Station Lid is **fully** open.
- b) To prepare the Injector, place it into the Reset Station. Push the Injector nose first into the Reset Station, ensuring the blue retractable core retracts, as shown in figure 1. Ensure the Injector is fully nested in the Reset Station.


figure 1

NOTE: Make sure the pin is in line with the Injector as shown in figure 2.

figure 2

- c) **FIRMLY** push on the lid of the Reset Station. This action is best performed from a standing position.
- d) Continue to push down firmly against the resistance until the lid closes completely, as shown in figure 3. The Lid will lift indicating the Injector is charged.


figure 3

NOTE: The Injector Hammer will be recessed when charged.

uncharged

charged

2) Fill Syringe

- a) Remove protective cap from the primary vial packaging and prepare the top following aseptic techniques.
- b) Aseptically remove the Filling/Vial Adapter from its packaging.
- c) Aseptically remove the Syringe from its packaging. Remove the blue protective Syringe Cap.
- d) Insert the Syringe into the top of the Filling Adapter. Push down and twist to the right (clockwise) to lock Syringe into place, as shown in figure 4. Then remove the protective sleeve from the Filling Adapter spike.
- e) If needed, remove the clear protective sleeve from the Filling Adapter, as shown in figure 5.
- f) Insert the Filling/Vial Adapter spike into the vial with a twisting motion until it is fully inserted, as shown in figure 6.
- g) Invert the assembly and pull down on the Syringe Plunger to draw fluid into the Syringe, as shown in figure 7. Remove air bubbles as required.
- h) Complete filling by aligning the red rubber ring with the Fill Line, as shown in figure 8.
- i) Snap off the end of the Plunger protruding from the Syringe, as shown in figure 9.

NOTE: The Filling/Vial Adapter is for single vial use. It may be used to continuously fill several Syringes from the same multi-dose vial for up to 48 hours. Up to five vial punctures with a new Filling/Vial Adapter can be performed on a multi-dose vial. Cap the Filling/Vial Adapter with a new unused Syringe between fillings.

NOTE: Prior to vial storage, remove and dispose of the Filling/Vial Adapter in accordance with local requirements.

figure 4

figure 5

figure 6

figure 7

figure 8

figure 9

3) Load Injector

NOTE: The Multi-function Button has two uses: 1) It locks the blue retractable core when it is fully forward while inserting a Syringe. 2) It activates the Injector to give an injection when the Syringe is pushed against the skin.

NOTE: To prevent inadvertent activation, the Multi-function Button should be depressed during steps 3a through 3c ensuring the Polished Core Ring is visible, as shown in figure 10.

figure 10

a) Ensure the Polished Core Ring on the blue retractable core is visible prior to depressing the Multi-function Button and loading a Syringe, as shown in figure 11.

figure 11

b) While depressing the Multi-function Button, insert the Syringe into the Injector until it “clicks” into place, as shown in figure 12.

figure 12

- c) Invert the assembly and remove the Syringe from the Filling/Vial Adapter by twisting to the left (counterclockwise), as shown in figure 13.

figure 13

NOTE: If injection is not to be given immediately, replace protective Syringe Cap.

4) Give Injection

NOTE: If needed, remove protective Syringe Cap.

- Remove thumb from the Multi-function Button, allowing the blue retractable core to move, and apply the Syringe perpendicular (90°) to the prepared injection site.
- Push the Syringe against the injection site until the blue retractable core stops moving.
- To inject, depress the Multi-function Button to activate the Injector. A “click” will be heard upon injection.
- Verify Syringe contents were delivered by inspecting that the Syringe Plunger has advanced completely to the Syringe tip.

Intramuscular Injection:

- Injecting straight into the tissue maximizes the depth of penetration, as shown in figure 14.

figure 14

Subcutaneous Injection:

- Pinch skin over the injection site to limit the depth of penetration, as shown in figure 15.

figure 15

5) Eject Syringe

- a) Eject the Syringe from the Injector by depressing the Syringe Ejection Button, as shown in figure 16.

figure 16

NOTE: The Syringe Ejection Button will move when pushing the Syringe against the skin. This is normal. Engage the Syringe Ejection Button only as shown to eject the Syringe. Push, do not slide or pull.

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