PharmaJet[®] Tropis[®]

EN – Needle-free Injection System
DE – Nadelloses Injektionssystem

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
CE
2797

Symbols Glossary

(Note: All referenced symbols are derived from **ISO 15223-1, Medical Devices - Symbols to be used with information to be supplied by the manufacturer.**)

Symbol	Title of Symbol	Description of Symbol	Ref. No.
	Manufacturer	Indicates the medical device manufacturer.	5.1.1
	Authorized representative in the European Community/ European Union	Indicates the authorized representative in the European Community/ European Union.	5.1.2
	Use-by date	Indicates the date after which the medical device is not to be used.	5.1.4
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	5.1.5
	Catalogue number	Indicates the manufacturer's catalogue number so the batch or lot can be identified.	5.1.6
	Date of manufacture / Country of manufacture	Indicates the date when the medical device was manufactured / to identify the country of manufacture of products.	5.1.3 / 5.1.11
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.	5.2.3
	Do not use if package is damaged	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.	5.2.8
	Single sterile barrier system	Indicates a single sterile barrier system.	5.2.11
	Single sterile barrier system with protective packaging outside	Indicates a single sterile barrier system with protective packaging outside.	5.2.14
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	5.3.7
	Do not re-use	Indicates a medical device that is intended for one single use only.	5.4.2
	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use.	5.4.3
	Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.	5.4.4
	Non-pyrogenic	Indicates a medical device that is non-pyrogenic.	5.6.3
	Medical device	Indicates the item is a medical device.	5.7.7
	Unique device identifier	Indicates a carrier that contains unique device identifier information.	5.7.10
	Do not touch	Indicates components that should not be touched.	N/A
	Prescription Only	Caution: US Federal law restricts this device to sale by or on the order of a physician.	N/A

PharmaJet[®]

Tropis[®]

English

System Components

Injector

Claimed Lifetime: 20,000 cycles

Syringe

STERILE EO

WARNING: Keep away from children of 3 years and younger.
Contains small parts.

Filling Adapter

STERILE EO

CAUTION: Federal law restricts this device to sale by or on the order of a physician.


Indications

Indications for Use: The PharmaJet Needle-free Injection System is intended to deliver various medications and vaccines intradermally by means of a narrow, high velocity fluid jet, which penetrates the skin and delivers the medicine or vaccine into the dermis. Healthcare providers who routinely administer injections may use the PharmaJet Needle-free Injection System. It may also be used by patients authorized by their physicians to self-inject, or have other individuals administer injections of prescribed medication.

Intended Users: The Tropis Needle-free Injection System (NFIS) is intended for use by healthcare providers or other authorized individuals and patients for administration of prescribed medication.

Intended Patient Population: It may be used for adults and children.

Intended Use Environment: The use environment is in home, in office, in clinic, and in the field.

NOTE: If the ambient temperature is outside of the specified range, DO NOT USE.

Use Specification

The PharmaJet Tropis NFIS can be used with medications intended for injection into the intradermal tissue. Types of medications that may be used with the Tropis NFIS include anti-infectives (e.g., vaccines), antineoplastic and immunomodulating agents, hormonal preparations, and other therapeutic products intended for 0.1mL intradermal injection. It is the responsibility of the prescribing physician to ensure medicines are used within their licensed indications and methods of administration.

Information to be Communicated to the Patient

- Some medications may produce a burning, aching, or stinging sensation. Check the medication information for a description of local reactions.
- Following the injection, a drop of blood may appear. Blood and/or fluid loss can be minimized by immediately applying pressure to the injection site for approximately 60 seconds with a cotton ball or gauze. Apply adhesive bandage if needed or desired.
- Patients receiving anticoagulants or persons with hemophilia or other coagulation disorders may have more postinjection bleeding and/or bruising.

Suitable Injection Sites

 **WARNING: DO NOT position the Syringe directly over a blood vessel. Although injection into the vessel is highly unlikely, injection over a vessel may produce an adverse reaction including, but not limited to, bleeding, bruising, swelling, or tenderness.**

Select a site that is free from irritation, with a low likelihood of infection, and where the skin is intact and free from abrasions, bruises, etc.

Needle-free injection sites are identical to those used for conventional needle and syringe; the deltoid, suprascapular, and lateral thigh areas are most often used as the delivery sites for vaccine. Other types of injectant may have specific sites for injection. Follow package directions from the manufacturer of injectant.

Safety and Clinical Performance / Clinical Benefits

The Tropis NFIS has been designed in accordance with applicable standards and regulations for safety and performance. The CE mark indicates PharmaJet products meet all relevant EU regulations. In addition, products meet WHO PQS requirements for quality and reliability under field immunization conditions with a safety profile ensuring no harm to patients, users or the environment. Adherence to PharmaJet's Quality Management System for manufacturing controls, design controls, and clinical evaluation assures favorable clinical benefit-to-risk ratio over the lifetime of the product.

Warranty Information

The Injector carries a replacement warranty against defects in materials and workmanship for two years from the date of purchase.

Contact Information

PharmaJet®

US Office:

PharmaJet, Inc.

400 Corporate Circle, Suite N

Golden, Colorado 80401

USA

Tel: +1 888 900 4321

Tel Outside the U.S.: +1 303 526 4278

EC REP

EU Authorized Representative:

Emergo Europe

Westervoortsedijk 60


6827 AT Arnhem

The Netherlands

Component Description


Injector

1. Activation Button
2. Inner Core Housing
3. Pin
4. Eject Button
5. Safety Button
6. Handles
7. Latch Button


Syringe

1. Syringe Cap
2. Syringe Barrel
3. Syringe Plunger
4. Fill Line

Filling Adapter

1. Syringe Interface
2. Spike

⚠ CAUTION: Use only PharmaJet components with the Tropis® Needle-free Injection System.

Warnings and Cautions

- To ensure delivered dose to target tissue, refer to **Workflow Section 4) Give Injection** for proper technique. For a partial or incomplete injection, refer to the medication or vaccine product labeling.
- In the event of changes in performance of the device, DO NOT USE, contact PharmaJet.
- Disposables (Syringe and Filling Adapters) components are single use ONLY. Reuse of components may lead to cross contamination.

⚠ CAUTION: Visually inspect sterile package for damage. DO NOT use if torn or damaged. DO NOT re-sterilize. Discard disposable components where the package is torn, damaged, or unintentionally opened prior to use.

Workflow

1) Prepare Injector

⚠ CAUTION: Repeated dry firing will damage Injector. Dry firing is activating the Injector with either an empty Syringe, partially-filled Syringe, or without a Syringe.

⚠ CAUTION: DO NOT open Injector Handles when device is charged.

- Verify if the Injector is in the charged or uncharged state by confirming the position of the Injector Pin as shown in figures 1a (uncharged) and 1b (charged). If the Injector is charged and the Handles are closed, proceed to **Section 2) Fill Syringe**. DO NOT open the Handles when the Injector is charged.

figure 1a

figure 1b

- When the Injector is in the uncharged state, as shown in figure 1a, and the Handles are closed, open the Handles to prepare the Injector. Press the rear Injector Latch Button to allow the Injector Handles to open freely, as shown in figure 2.


figure 2

NOTE: If the Injector is in the charged state, as shown in figure 1b, and the Handles are partially open, as shown in figure 3, close the Handles and proceed to **Section 2) Fill Syringe**. DO NOT touch internal components or attempt to fully open the Handles when the Injector is charged.


figure 3

c) With the Injector Handles fully open, as shown in figure 4, squeeze the Handles until they close and are locked in place, as shown in figure 5. Once the Injector Handles are closed, ensure the Injector Pin is in the charged state, as shown in figure 1b.

figure 4

figure 5

NOTE: If the Injector continuously fails to reset, stop use of the device. Refer to the contact information found at the end of this document and contact a PharmaJet representative.

2) Fill Syringe

WARNING: The tip of the Filling Adapter spike is sharp. Use caution when handling.

WARNING: Syringes auto-disable during use. DO NOT reuse a Syringe, as it may result in contamination or device malfunction.

⚠ WARNING: Vial integrity or sterility is not ensured when used with the Filling Adapter. DO NOT reuse the Filling Adapter in another vial. Doing so may contaminate contents of new vial.

⚠ CAUTION: The protective Syringe Cap is designed to protect against touch contamination only, and does not provide a barrier against other forms of contamination.

⚠ CAUTION: DO NOT attempt to use a partially filled Syringe; damage may occur.

- a) Remove protective cap from the primary vial packaging and prepare the top of the vial following aseptic techniques.
- b) Aseptically remove the Filling Adapter from its packaging.
- c) Aseptically remove the Syringe from its packaging. Remove Protective Cap.
- d) Insert the Syringe into the top of the Filling Adapter. Push down and twist to the right (clockwise) to lock Syringe into place, as shown in figure 6. Then remove the protective sleeve from the Filling Adapter spike, if present.
- e) Insert the Filling Adapter spike into the vial with a twisting motion until it is fully inserted, as shown in figure 7.

NOTE: Adding air to a vial can minimize vacuum and assist filling.

- f) Invert the assembly.
- g) Pull down Plunger to draw fluid into Syringe. Remove bubbles as required. Complete filling by aligning red rubber ring with the fill line, as shown in figure 8.
- h) Snap off end of Plunger protruding from Syringe, as shown in figure 9.

NOTE: The Filling Adapter is for use on one vial only. Once the vial is punctured, the Filling/Vial Adapter must remain in place for the life of the vial. A new (unused) Syringe must be placed on the Filling Adapter at all times; immediately after or between injections and when the vial is in storage during and after immunization sessions. If at any time and for any reason the Filling Adapter is removed, the vial must be discarded.

figure 6

figure 7

figure 8

figure 9

3) Load Injector

⚠ WARNING: DO NOT direct the Syringe toward the face or nearby skin when loading the Syringe into an Injector, preparing to give an injection, or ejecting a Syringe after use. Injury may occur.

⚠ CAUTION: Partial or incorrect insertion of a Syringe into an Injector may damage the Syringe and/or Injector.

⚠ CAUTION: DO NOT load a Syringe if the Inner Core Housing is retracted. This may cause inadvertent activation.

- a) Insert the Syringe and Vial Assembly into the Inner Core Housing of the charged Injector until it “clicks” into place, as shown in figure 10.


figure 10

- b) Invert the assembly and remove the Syringe from the Filling Adapter by twisting to the left (counterclockwise).

NOTE: If injection is not to be given immediately, replace protective Syringe Cap.

4) Give Injection

⚠ WARNING: Improper injection site location or technique may result in injury.

⚠ CAUTION: If the Injector fails to reset or if the Inner Core Housing does not move freely, DO NOT USE. Continued use may result in an incomplete injection or failure to inject.

⚠ CAUTION: Patients receiving anticoagulants or persons with hemophilia or other coagulation disorders may have more post-injection bleeding and/or bruising.

⚠ CAUTION: DO NOT deliver injection with the Injector at an angle other than 90 degrees to the patient, as shown in figure 12.

NOTE: If needed, remove protective Syringe Cap.

a) With forefinger, depress and hold the yellow Injector Safety Button, as shown in figure 11. This allows the Inner Core Housing to retract freely.


figure 11

b) While holding the Safety Button, press the Syringe perpendicular to the injection site, as shown in figure 12. Once the Inner Core Housing is retracted, it is not necessary to hold the Safety Button as long as you keep pressure on the patient and the Inner Core Housing remains retracted.

figure 12

c) To inject, depress Activation Button with thumb. A “click” will be heard upon injection.

d) Verify Syringe contents were delivered by observing that the Plunger has advanced completely to the Syringe tip.

NOTE: If the Inner Core Housing returns to the original extended position and becomes locked, simply press the Safety Button again (See Figure 11) and start injection operation from Section 4.

5) Eject the Syringe

a) Remove Syringe from the Injector by pressing the Eject button, as shown in figure 13.


figure 13

- b) Dispose of Syringe in accordance with local requirements.
- c) Visually inspect the Injector for contamination of foreign particles and soil. If device is visibly contaminated, perform cleaning and disinfection according to Cleaning and Disinfection Instructions.
- d) For storage, keep the Injector in the uncharged and closed configuration.

Cleaning, Disinfecting, Storage, and Disposal

CAUTION: DO NOT autoclave or sterilize the Injector or disposables (Syringe and Filling Adapter) components.

INSPECTION: Visually inspect the exposed surfaces of the opened Injector for foreign contaminants and soil. If foreign contaminants or soil are present, proceed to cleaning and disinfection. Do not disassemble the device for inspection or for cleaning and disinfection.

CLEANING: Rinse the Injector with potable tap water for at least 2 minutes, rubbing and wiping large debris off device. Wipe device with 70% Isopropyl Alcohol to remove visible soil. Dry the Injector with a clean, absorbent cloth. Visually inspect Injector after cleaning. If debris remains, repeat procedure.


DISINFECTING: Spray the Injector with 70% Isopropyl Alcohol and wipe exposed surfaces for at least 1 minute. Maintain surface wetness with 70% Isopropyl Alcohol during disinfection process.

NOTE: *Injectors should be periodically cleaned and disinfected to maintain device life.*

STORAGE: It is recommended to store the Injector in the uncharged and closed configuration, in the protective case.

DISPOSAL: Dispose of used devices in accordance with local requirements.

**Storage and Transport
Temperature:**

Reporting of Serious Incidents

The user or the patient shall report any serious incidents associated with the PharmaJet Tropis NFIS to PharmaJet and to appropriate regulatory bodies.

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