

Modulating immune responses to personalized cancer vaccines using PharmaJet Needle-free Injection Systems


Carmen Ledesma-Feliciano, Kira Elma, Shaily Garg, Nathalie Landry - PharmaJet (Golden, CO)

PharmaJet® is transforming vaccine delivery through innovative Needle-free Injection Systems. Our enabling technology can enhance immune response and provides safety, efficacy, efficiency, and patient comfort using intradermal (ID), intramuscular (IM), or subcutaneous (SC) applications.^{a, b, c}

PharmaJet®



Tropis® ID
provides intradermal injections



Stratis® IM/SC
provides intramuscular or subcutaneous injections

Recent Strategies for Personalized Cancer Vaccines


Personalized vaccine strategies for cancer treatment include nucleic acid platforms with antigen presenting cell (APC) targeting to boost T cell activation and the addition of neoantigenic epitopes. Combined with checkpoint inhibitors (CPIs), therapeutic cancer vaccines significantly enhance anti-tumor and clinical responses. PharmaJet devices can improve DNA vaccine delivery^d and have been successfully adopted into multiple vaccine development programs while demonstrating both immunogenicity and clinical efficacy for novel personalized therapeutic vaccines as shown below.



VB10.NEO DNA vaccine Phase 1/2a Trial


INTRODUCTION: VB10.NEO was administered to patients with locally advanced or metastatic solid cancers in combination with CPI and/or other anti-cancer therapies at investigator's discretion.

Each individual VB10.NEO vaccine contains up to 20 patient-specific neoantigens and is designed to target APCs (CCL3L1) using Nykode Therapeutics' modular Vaccibody™ vaccine platform.




VB10.NEO administration with Stratis:

T cell responses (n=7, 22 weeks)¹



Durability of De Novo Vaccine Responses (n=8, 100 weeks)²



CONCLUSION: Assessment of neoantigen-specific T cell reactivity demonstrated VB10.NEO induced **broad and long-lasting T cell responses** almost a year after the last vaccination, and the majority of tested neoantigens activated polyfunctional CD8+ T cells.


¹Krauss, et al (2023). Individualized APC targeting VB10.NEO cancer vaccines induce broad neopeptide-specific CD8 T cell responses in patients with advanced or metastatic solid tumors: interim results from a phase 1/2a trial. American Association for Cancer Research (AACR) 2023 poster. ²Weinert, et al (2025). DNA-encoded individualized neoantigen vaccines elicit durable and functional immune responses. Association for Cancer Immunotherapy CIMT 2025 poster.




EVX-02 & EVX-03 DNA vaccines Phase 1/2 Trial

INTRODUCTION: EVX-02 was administered to patients that had a complete resection of Stage IIIB/IIIC/IIID or Stage IV melanoma and were at a high risk of recurrence, in combination with CPI (nivolumab).


Each EVX-02 vaccine encoded up to 13 patient-specific neoantigens. Study compared Stratis delivery to poloxamer delivery with needle and syringe (NS). EVX-03 additionally encodes APC-targeting and dimerization units.



EVX-02 administration with Stratis (n=5):



- Induced neoantigen-specific T cell responses in all patients and CD4+ and CD8+ T cells contributed to these responses.
- Immune responses trended higher in response magnitude with Stratis (cohort B) compared to poloxamer delivery with NS (cohort A).



EVX-03 delivered by PharmaJet in a mouse pre-clinical model induced a strong antitumor effect and neopeptide-specific T cell response more potent than 1st generation delivery with poloxamer.

CONCLUSION: EVX-02 was immunogenic and efficacious when delivered with Stratis, showing an improvement in T cell induction compared to NS. These patients were relapse-free at their last assessment (1 year after treatment start). The 2nd generation EVX-03 vaccine induced a more potent response in a mouse pre-clinical model.


Kleine-Kohlbrecher, et al (2023). AI-designed personalized neoantigen vaccine, EVX-02, induces robust T-cell responses in melanoma patients. Society for Immunotherapy of Cancer (SITC) 2023 poster.



SCIB1/ iSCIB1+ DNA vaccines Phase 2 Trial


INTRODUCTION: The SCOPE trial is evaluating SCIB1 and next-generation iSCIB1+ therapeutic vaccines delivered with Stratis IM, combined with checkpoint inhibitors in patients with advanced melanoma. A new trial arm is evaluating Tropis ID delivery of iSCIB1+.

SCIB1 (Cohort 1): ImmunoBody® DNA vaccine, TRP-2, gp100; limited to specific human leukocyte antigen (HLA) types.
iSCIB1+ (Cohort 3): AvidiMab®-modified, includes additional melanoma epitopes, broadening HLA coverage.



SCIB1 and iSCIB1+ administration with Stratis^{1,2}:

Cohort 1 & 3 RECIST 1.1 Best Overall Response (BOR): Waterfall Plot



T cell Responses - ELISpot Assays (n=33)

CLINICAL RESPONSE	NUMBER OF PATIENTS	POSTIVE T CELL RESPONSE	% POSITIVE
CR/PR	19	15	79%
SD	7	5	71%
PD	7	4	57%


Summary

- Both Cohort 1 and Cohort 3 showed impressive efficacy in the proposed development population.
- CD8 T cell responses are important for clinical responses.
- The safety profile of SCIB1& iSCIB1+ is benign with no potentiation of CPI toxicities (data not shown).
- iSCIB1+ selected for further development with PharmaJet delivery: shows eqi-potency and safety compared to SCIB1 and is efficacious in a wider population.

CONCLUSION: SCIB1 and iSCIB1+ DNA vaccines are immunogenic and lead to significant melanoma lesion size regression in a majority of treated patients.

¹Scancell (2025, July 22). Active Immunotherapy for a Cancer-Free Future: Immunobody® iSCIB1+ strongly improved outcomes in Late-Stage Melanoma SCOPE study results. [Company presentation]. ²Scancell (2025, September). Active Immunotherapy for a Cancer-Free Future [Company presentation].

Best T cell Response (Cohort 3)



- All six iSCIB1+ epitopes generated CD8+ T cell responses.
- 72% of patients responded to both TRP-2 & gp100, reducing risk of immune escape.
- 19/31 (61%) HLA-matched patients made a T cell response
- 15/19 (79%) amongst clinical responders.
- 10/12 (83%) with a CD8+ T cell response were clinical responders.

Conclusions

- PharmaJet devices enable DNA personalized cancer vaccine delivery: Vaccinations are safe and well tolerated, induce neoantigen-specific T cell responses, and lead to favorable clinical outcomes.
- PharmaJet devices can be easily incorporated into novel personalized therapeutic strategies to treat various types of cancer.

References: ^aAlberer, M. et al. (2017). Safety and immunogenicity of a mRNA rabies vaccine in healthy adults: an open-label, non-randomised, prospective, first-in-human phase 1 clinical trial. *Lancet* 390, 1511-1520. ^bMohan, D. et al. (2025). Evaluating the Impact of Needle-Free Delivery of Inactivated Polio Vaccine on Nigeria's Routine Immunization Program: An Implementation Hybrid Trial. *Vaccines (Basel)* 13. ^cBavdekar, A. et al. (2018). Immunogenicity and safety of measles-mumps-rubella vaccine delivered by disposable-syringe jet injector in India: A randomized, parallel group, non-inferiority trial. *Vaccine* 36, 1220-1226. ^dLedesma-Feliciano C, et al (2023). Improved DNA Vaccine Delivery with Needle-Free Injection Systems. *Vaccines*. 11(2):280. Poster presented at PCVS on November 19, 2025. Doc #60-10727-007G